Liebig-Dispenser

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Diffuseur-Liebig

Mode d'emploi

Diffusore-Liebig

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Liebig-Dispenser

Instructions

Liebig verdamper

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
empire – applicare – controllare

EN For the treatment of varroosis of honey bees in combination with formic acid in concentrations of 60% or 85%. Best suited for use together with the veterinary medicinal products Formivar 60% or Formivar 85%.

Components

- ① Base plate with reservoir and fixing prongs for a secure fixation of the bottle.
- ② Perforated paper-wick for easy and exact dosage.
- ③ Dropper-insert with calibrated outflow for constant dosage.
- ④ Graduated low-profile bottle for easy filling and quick control during treatment.

Mode of action

The formic acid is delivered slowly and steadily through the calibrated opening of the dropper into the reservoir of the base plate. By capillary action, the acid is distributed evenly over the entire surface of the paper-wick. The intensity of the evaporation is determined by the size of the paper-wick.

The paper can be folded along the fine perforated lines and the excess paper can be accurately and neatly removed. As long as there is formic acid in the bottle, the evaporation will be constant. The treatment is finished as soon as the paper-wick is dry again.

Sizing of paper-wick

For the different concentrations of formic acid, follow the instructions of the paper-wick, green for 60% and red for 85%. On the paper-wick you will find the information on the recommended dosage for the first and second treatment which depend on the hive system and the temperature.

Be aware of the expected maximum daytime temperature during the first days of treatment. The size of the paper-wick is reduced by folding and separating the paper along the perforations as demonstrated in the illustrations by the dark shaded areas. The amount of formic acid evaporation is regulated by the size of the paper-wick. A variety of external factors such as the type and size of the hive, outside temperature and sunlight can influence the amount of evaporation. The general instructions of the paper-wick should serve as a guide. For example, for Formivar 60%, 20–50 ml of formic acid can evaporate per day depending on the hive.

Additional information

The formic acid treatment with the Liebig-Dispenser has proven to be safe for the bees for single and double brood chamber hives as well as for Swiss hives. Young brood and emerging bees are very susceptible to formic acid which is why a formic acid treatment can always lead to some brood losses or brood stop. These losses have no negative influence on the later development of the colony.

When treating with Formivar 85% and an outside temperature of 25°C or above, the colony may expand and take up more space to improve the control of the formic acid concentration.

In some cases, bearding up of the bees in front of the hive entrance may occur. It is therefore important that there is sufficient space for the expansion of strong colonies. We recommend that you take additional measures so that continuous ventilation is guaranteed. For example, open the varroa sliding tray one third or switch to formic acid 60%.

Cleaning

Remove propolis with HalaApi 898. Clean the base plate from wax residues with Apisoda in warm water. Do not clean the plate in the dishwasher as it will warp at temperatures above 40°C.

Treatment instructions


1. Remove the dropper-insert.
2. Fill the bottle up to the marking with the recommended amount.
3. Insert the dropper. Close the bottle for transport with the childproof cap and store at ambient temperature.
4. Size the paper-wick according to the type of hive and the weather forecast.
5. Position the base plate over the center of the brood frames and place the adjusted paper-wick. The base plate should lie as stable and horizontal as possible so that the reservoir does not overflow.

- Place the bottle upside down on the fixing prongs without applying pressure. The formic acid then drips into the flat reservoir and is distributed throughout the paper. **Caution:** Only use acid that is of ambient temperature as acid that is colder will expand and potentially overflow.
- Check that the formic acid is flowing out properly by slightly tilting the bottle.
- For full size colonies, open the hive entrance completely to improve ventilation. For small colonies such as nuclei and splits, leave the hive entrance as it is to reduce the chance of robbing.
- Check the evaporated amount after one day at the earliest and adjust the size of the paper-wick if necessary; we recommend another check after 3 days at the latest. As soon as the bottle is empty and the paper is dry, remove the Liebig-Dispenser from the hive.

Important

- Close the floorboard at least 3 days before treatment so that the colony can adjust to the reduced ventilation.
- Do not reduce or rearrange the frames in the brood chamber before treatment. Allow sufficient space for the colony to expand for ventilation.
- The dispenser should be placed as horizontally as possible.
- Start the treatment preferably early in the morning if high temperatures are expected.
- To reduce bee stress, do not start treatment shortly before thunderstorms or heavy rain.
- For unaltered honey, start treatment only after the honey harvest.
- Do not treat and feed simultaneously. The formic acid binds with fresh feed and there is also a risk of robbery.
- 1st treatment: end of July to mid-August before feeding.
- 2nd treatment: after feeding in September.
- We recommend a new paper-wick for each treatment.

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Andermatt BioVet AG